

Epiretinal Membrane and Instrument Segmentation in Retinal Surgery Video

Saurav Verma


College of Information Science, University of Arizona
Advisor: Dr. Eungjoo Lee

Vision, Systems & Intelligence (VSI) Lab — University of Arizona

PURPOSE


- ERM peeling depends on accurately identifying thin membranes and instrument-tissue interaction.
- Ophthalmic surgical video remains underrepresented in pixel-level segmentation research.
- This project develops a pilot dataset and evaluates modern deep-learning models for ERM and tool segmentation.

CLINICAL CONTEXT



Intraoperative view showing forceps initiating an ERM peel.

ANNOTATION EXAMPLES



Annotated Overlay + Greyscale Mask


ANNOTATION SUMMARY

- 300 frames (1 fps) from 7 ERM surgery videos.
- Annotations done in *Supervisely* framework.
- Classes: ERM, Retinal Forceps, Light Tool.
- Surgery phase labels (Flap initiation, Dye, Peeling, Completion) collected for future modeling.
- Includes challenging cases: faint ERM, tool overlap, illumination artifacts.

DATASET & PREPROCESSING


- Frame extraction, deduplication, blur filtering.
- Region-of-Interest cropping to isolate surgical field.
- Resize: 512×512, 768×768.
- Surgery-level train/val/test split.

AUGMENTATION PIPELINE



Training augmentations: rotation, flips, color jitter, scaling
(Top-left: original frame/mask; other panels show a few augmented variants)

END-TO-END WORKFLOW




QUANTITATIVE RESULTS

Model	ERM	Forceps	Light Tool (Endoillumination)
DeepLabV3+	0.77	0.83	0.76
Unet++	0.70	0.80	0.72

Dice Scores for ERM & Instruments
(This is an ongoing research project; detailed internal evaluations are not disclosed)

QUALITATIVE RESULTS



Top Row: Successful segmentation of ERM + instruments.
Bottom Row: Challenging case : artifacts glare, faint ERM.
Column : Left: Original, Middle: GT, Right: Prediction

KEY INSIGHTS

- ERM & tool segmentation is feasible with limited data.
- 512 px resolution provides strong segmentation and contour accuracy despite extreme foreground sparsity (<1% pixels).
- DeepLabv3+ provides the most stable performance in this pilot dataset.

CONCLUSION

DeepLabv3+ at 512 px showed the most consistent performance, demonstrating that ERM and instrument segmentation from intraoperative video is technically feasible even with limited curated data.

FUTURE EXTENSIONS

- Expand dataset with more surgeries.
- Explore transformer-based architectures.
- Add surgical-phase recognition in a multi-task pipeline.